

APERTURE	MESH NO.	Number	APERTURE	MESH NO.	Number
2.80 mm	7	31	5 Micron	1400	1
3.35 mm	6	32	10 Micron	1000	2
4.00 mm	5	33	20 Micron	635	3
4.75 mm	4	34	25 Micron	500	4
5.60 mm	3,1/2	35	37 Micron	400	5
6.35 mm	1/4"	36	45 Micron	325	6
8.00 mm	5/16"	37	53 Micron	270	7
9.51 mm	3/8"	38	63 Micron	230	8
11.20 mm	7/16"	39	75 Micron	200	9
12.70 mm	1/2"	40	90 Micron	170	10
16.00 mm	5/8"	41	105 Micron	150	11
19.00 mm	3/4"	42	106 Micron	140	12
22.60 mm	7/8"	43	125 Micron	120	13
25.40 mm	1"	44	150 Micron	100	14
31.50 mm	1,1/4"	45	180 Micron	80	15
38.10 mm	1,1/2"	46	212 Micron	70	16
44.50 mm	1,3/4"	47	250 Micron	60	17
50.80 mm	2"	48	300 Micron	50	18
63.50 mm	2,1/2"	49	355 Micron	45	19
76.20 mm	3"	50	420 Micron	40	20
89.00 mm	3,1/2"	51	500 Micron	35	21
101.60 mm	4"	52	600 Micron	30	22
			710 Micron	25	23
—	Washing sieve 200	53	840 Micron	20	24
—	2×20 sieve	54	1.00 mm	18	25
—	1×20 sieve	55	1.20 mm	16	26
—	Pan of sieve	56	1.40 mm	14	27
—	Cover of sieve	57	1.70 mm	12	28
			2.00 mm	10	29
			2.36 mm	8	30